



VALVES FOR HIGHEST REQUIREMENTS



A vision becomes reality

AVAX is the trend-setting and unrivalled split butterfly valve technology from Andocksysteme. Completely made of medical-grade plastic materials it represents the origin for a new generation of single-use containment interfaces.

Strictly designed as a congruent passive to passive system, this unique valve enables its usability in any desirable application due to its interoperability.



Ready to connect



Docked and locked



- Meets OEB 5 requirements (STTW ≤ 1 mcg/m³)
- Intuitive and secure handling
- Tri-Clamp seal integrated (4" BS 4825-3)
- Hybrid operation with stain less steel AVAX optional
- Sterile execution available



Key technology of AVAX is the implementation of an established locking mechanism this technical evolution it is the first time that housings and discs are locked to each other individually but by a mutual locking motion.



Locking features for housing and disc

As a result of two component injection moulding manufacturing process, fully hygienic is guaranteed. The seal material is incorporated by chemical bonding forces wich gives seamless results. All system components comply with the major pharmaceutical requiréments like FDA, USP and EP.





via AVAX.

Safe filter change



In transfer mode. Integrated seals



Protective cap and originality under seal



Multifunctional lever used as a handle



Are you still cleaning or being productive



Contained sampling

Andocksysteme G. Untch GmbH Schweighofstr. 3 79410 Badenweiler Tel. +49 7632 82 3 82 0 Fax +49 7632 82 3 82 28 info@andocksysteme.de

Sales Representative